



A. D. Darbal^{1,4}, R. D. Narayan¹, C. Vartuli², G. Lian², R. Graham³, F. Shaapur³, S. Nicolopoulos⁴ and J. K. Weiss¹

¹: AppFive LLC, Tempe, Arizona, USA. ²: Texas Instruments, Dallas, Texas, USA ³: Nanotem Inc, Scottsdale, Arizona, USA ⁴: NanoMEGAS SPRL, Brussels, Belgium

PRECESSION ELECTRON DIFFRACTION



Portillo, J., et al. (2010). *Materials Science Forum* **23**, 1-7




Diffraction patterns from two points 120 nm apart from Si/SiGe multilayered specimen

- No particular beam is strongly diffracted – reducing strong dynamical effects
- Insensitive to small thickness and orientation changes
- Number of spots increases – better sampling of higher order spots

STRAIN MEASUREMENT ANALYSIS

- Diffraction patterns from strained region are matched against a reference pattern.
- Reference pattern from unstrained region.
- Correlation distance used as the metric for fitting reference to strained patterns.
- Results include strain in x and y-directions and shear (not shown).
– Relative to x-direction specified by user.



INTRODUCTION

- Measurement of strain with high spatial resolution and high precision in semiconductor devices is critical to monitor the designed and unintended strain distributions.
- Use of spot diffraction patterns with nanobeam illumination gives higher spatial resolution than other TEM techniques[1].
– Experiment is relatively simple.
- Technique is made possible by beam precession
– Improves quality of diffraction patterns..

[1] D Cooper et al., *Journal of Physics: Conference Series* **326** 012025 (2011).

LIMITATIONS OF EXISTING TEM STRAIN MEASUREMENT METHODS

Technique	Advantages	Limitations
Convergent beam electron diffraction	High spatial, strain sensitivity	<ul style="list-style-type: none"> • Needs the sample to be relatively thick (>150 nm) • Sample needs to be oriented away from a low index axis. • Very sensitive to strain relaxation
Dark field holography	High spatial resolution (5 nm), large field of view (1 μm×1 μm)	<ul style="list-style-type: none"> • Requires unstrained reference with identical crystallographic orientation area close to strained region
High resolution imaging	High spatial resolution (< 1 nm)	<ul style="list-style-type: none"> • Limited field of view (100 nm² × 100 nm²) • Stringent requirements on specimen quality


CONVENTIONAL NANOBEBAM DIFFRACTION

- Acquire spot diffraction patterns from strained and unstrained regions using a quasi-parallel nanoprobe (<5 nm)
- Use measured shift in spot positions to calculate strain
- Experiment is relatively straightforward
- Limitations
 - Presence of strong dynamical effects lead to rapid changes in spot intensities with small thickness and orientation changes
 - Strong dependence of spot intensities on changes in local thickness and orientation makes automated analysis challenging
– Requires manual intervention in identifying spot positions
 - Inadequate sampling of higher order reflections limits the accuracy

RESULTS


- Data acquired with Zeiss Libra L200 TEM.
 - Field Emission Gun (FEG)
 - Scanning TEM (STEM) mode.
 - NanoMEGAS DigiSTAR unit for precession and descanning of beam.
- Positive percentage strain values correspond to tensile strain, negative values compressive.
- Precision of strain measurement is 0.02% in profile below.
– Precision can depend on characteristics of data, such as resolution of spot patterns.

1. Strain profile of an Si/SiGe layer.




- Strain in x-direction is near zero, indicating a coherent interface.

2. Strain maps from a AlGaIn/GaN HEMT. (provided by Cree)



- For such devices, tensile strain is expected to be asymmetric on different sides of the gate, and this is seen here.
- Compressive strain perpendicular to the interface (x-direction) of ~1% is seen in the AlGaIn region.

3. Strain maps from the Si region of a pMOS device.



- x and y-directions aligned with [220] and [002] directions in Si.
- Localized biaxial tensile strain close to contact edges.